

Cambridge International Examinations Cambridge International Advanced Subsidiary and Advanced Level

BIOLOGY

Advanced Practical Skills 1

9700/33 October/November 2014

CONFIDENTIAL INSTRUCTIONS

Great care should be taken to ensure that any confidential information given, including the identity of material on microscope slides where appropriate, does not reach the candidates either directly or indirectly.



If you have any problems or queries regarding these Instructions, please contact CIEby e-mail:info@cie.org.ukby phone:+44 1223 553554by fax:+44 1223 553558stating the Centre number, the nature of the query and the syllabus number quoted above.

This document consists of 8 printed pages.

DC (SJF/CGW) 70249/3 © UCLES 2014 **CAMBRIDGE** International Examinations

[Turn over

Instructions for preparing apparatus

These instructions give details of the apparatus required by each candidate for each experiment in this paper. A summary of the questions that will be presented to the candidates is included, where appropriate, to allow the biology teacher to test the apparatus appropriately.

No access to the question paper is permitted in advance of the examination.

If a candidate breaks any of the apparatus, or loses any of the material supplied, the matter should be rectified and a note made in the Supervisor's Report.

Candidates must be provided with a microscope with:

- Eyepiece lens, $\times 10$ (equal to 16 mm or $\frac{2}{3}$ ")
- Low-power objective lens, $\times 10$ (equal to 16 mm or $\frac{2}{3}$ ")
- High-power objective lens, $\times 40$ (equal to 4 mm or $\frac{1}{6}$ ")
- Eyepiece graticule fitted within the eyepiece and visible in focus at the same time as the specimen.

To avoid confusion, Cambridge request that only the lenses specified above are fitted in the microscopes to be used in the examination. Any lenses which are **not** $\times 10$ or $\times 40$ should be removed or replaced.

Each candidate must have sole, uninterrupted, use of the microscope for at least 55 minutes.


Supervisors are advised to remind candidates that **all** substances in the examination should be treated with caution. Pipette fillers and safety goggles should be used where necessary.

In accordance with the COSHH (Control of Substances Hazardous to Health) Regulations, operative in the UK, a hazard appraisal of the examination has been carried out.

The following codes are used where relevant.

- **C** = corrosive substance
- H = harmful or irritating substance
- **T** = toxic substance

- **F** = highly flammable substance
- **O** = oxidising substance
- \mathbf{N} = harmful to environment



Centres are reminded that they are **not** permitted to open the question paper envelopes before the examination. Centres should also refer to the Handbook for Centres.

If there are any difficulties with any aspect of setting up this practical examination that the Centre is not able to resolve, it is essential for Centres to contact the Product Manager as soon as possible by **e-mail** to <u>info@cie.org.uk</u>, by **fax** to +44 1223 553558 or by **phone** to +44 1223 553554.

Confidential Instructions

No access to the Question Paper is permitted in advance of the examination.

Each candidate will require:

For both Questions:

• mm ruler.

Question 1

- Solutions provided to the candidates should be supplied in a suitable beaker, or container, for removal of the solution using a syringe. More of the solutions should be available if requested by candidates.
- All solutions should be provided to candidates at **room temperature**.
- Fresh test-tubes and syringes are needed for each candidate.
- Fresh **S1**, **S2**, **P** and **W** are needed for each candidate.
- All solutions should be disposed of according to local safety regulations.

Summary of solutions:

labelled	contents	hazard	volume /cm ³
S1	0.10 mol dm ⁻³ sucrose solution	none	200
S2	1.00 mol dm ⁻³ sucrose solution	none	200
Р	0.50 mol dm ⁻³ sucrose solution coloured with methylene blue	none	20
W	distilled water	none	200

It is advisable to wear safety glasses/goggles when handling chemicals.

Preparation of solutions:

Sucrose solutions may be prepared the day before but should be kept in covered containers in a refrigerator. However, they should be supplied to candidates at **room temperature** for the examination.

(i) S1, at least 200 cm^3 of 0.10 mol dm^{-3} sucrose solution, in a beaker or container labelled S1.

This is prepared by dissolving 6.8g of sucrose in 100 cm³ of distilled water and making up to 200 cm³ with distilled water.

This is sufficient for 1 candidate.

(ii) S2, at least 200 cm^3 of 1.00 mol dm^{-3} sucrose solution, in a beaker or container labelled S2.

This is prepared by dissolving 68g of sucrose in 100 cm³ of distilled water and making up to 200 cm³ with distilled water.

This is sufficient for 1 candidate.

(iii) **P**, at least 20 cm³ of 0.50 mol dm⁻³ sucrose solution coloured with methylene blue, in a beaker or container labelled **P**.

This is prepared by dissolving 34 g of sucrose in 100 cm³ of distilled water and making up to 200 cm³ with distilled water.

Put 1 cm³ of 1% methylene blue solution into this 200 cm^3 of **P**, mixing well, so that it is coloured dark blue.

This is sufficient for 10 candidates.

1% methylene blue solution is prepared by dissolving 1g of methylene blue and 6g sodium chloride in 80 cm³ of distilled water and making up to 100 cm³ with distilled water.

(iv) W, at least 200 cm^3 of distilled water in a beaker or container, labelled W.

This is sufficient for 1 candidate.

Apparatus for each group of candidates should be clean. Syringe needles are **not** required and must **not** be given to candidates.

Apparatus for each candidate		1
10 cm ³ syringe with the means to wash it out		
$5{\rm cm}^3$ syringe with the means to wash it out – when put into the test-tube the syringe should be supported on the rim of the test-tube		
Glass rod	1	
Beaker or container(s), (about 400 cm ³), with about 300 cm ³ tap water for pouring into test-tubes, labelled For washing	2	
Container(s), to hold about 400 cm ³ , labelled For waste		
White card or white paper, approximately 100 mm by 100 mm		
Paper towels	8	
Test-tube – to hold a volume more than $40 \mathrm{cm^3}$ and less than $50 \mathrm{cm^3}$		
Beakers or containers to hold at least 120 cm ³		
Test-tube rack or container to hold the large test-tube		
Glass marker pen		
Safety goggles/glasses	1	

During the examination, the Supervisor (**not** the Invigilator) should, **out of the sight of the candidates**, carry out **Question 1** using the same solutions as the candidates. These results should be written in the Supervisor's Report, not on a spare question paper.

The Supervisor's Report and the candidates' seating plan should be enclosed with the candidates' scripts.

Please ensure that, if the scripts are in several packets, a copy of the Supervisor's Report and the candidates' seating plan are enclosed with each packet of scripts.

The Invigilator should **not** carry out **Question 1**.

https://xtremepape.rs/

Question 2

On receipt of the slides, please check that they are labelled **K1** and that none of the slides are broken. The material is **confidential** (so **must not** be disclosed to candidates) and the slides should **not** be viewed in advance of the examination.

The number of slides supplied by Cambridge will be equal to half the candidate entry.

Therefore, half the candidates should start on **Question 2** and the other half should start on **Question 1**.

For each candidate:

- the microscope must be set up on low power
- the slide must **not** be on the stage of the microscope.

Each candidate must have sole, uninterrupted use of the microscope for 55 minutes.

- (i) Slide K1 (supplied by Cambridge).
- (ii) Microscope with:
 - Low-power objective lens, ×10 (equal to 16 mm or ²/₃)
 - High-power objective lens, $\times 40$ (equal to 4 mm or $\frac{1}{6}$)
 - Eyepiece lens, $\times 10$ (equal to 16 mm or $\frac{2}{3}$ ")
 - Eyepiece graticule fitted within the eyepiece and visible in focus at the same time as the specimen.

To avoid confusion, Cambridge request that only the lenses specified above are fitted in the microscopes to be used in the examination. Any lenses which are **not** \times 10 or \times 40 should be removed or replaced.

MATERIALS TO BE SUPPLIED BY CAMBRIDGE

- (i) Question papers.
- (ii) Slide K1.

RETURN OF EXAMINATION MATERIALS TO CAMBRIDGE

Immediately after the examination the microscope slides must be:

- returned to Cambridge in the containers in which they were received, using the self-adhesive label. They must **not** be included in the packet of scripts.
- or
- purchased using the order form enclosed with the slides, which should be completed and returned to Cambridge. The order form must **not** be included in the packet of scripts.

Slides and boxes will be charged at the rate of £3 per slide and £1 per box.

If the items are not returned or purchased by the deadline stated on the order form they will be charged at £3.50 per slide plus £1 per box.

SUPERVISOR'S REPORT and SEATING PLAN

The teacher responsible for the examination is asked to fill in the Supervisor's Report in these Confidential Instructions. For Centres where more than one script packet is used, there must be a copy of the completed Supervisor's Report and the candidates' seating plan in each script packet.

These Supervisors' Reports are essential in order to allow the Examiners to assess all candidates as fairly as possible and should always be completed by every Centre.

https://xtremepape.rs/

Permission to reproduce items where third-party owned material protected by copyright is included has been sought and cleared where possible. Every reasonable effort has been made by the publisher (UCLES) to trace copyright holders, but if any items requiring clearance have unwittingly been included, the publisher will be pleased to make amends at the earliest possible opportunity.

Cambridge International Examinations is part of the Cambridge Assessment Group. Cambridge Assessment is the brand name of University of Cambridge Local Examinations Syndicate (UCLES), which is itself a department of the University of Cambridge.

This form should be completed and sent to the Examiner with the scripts.

7

SUPERVISOR'S REPORT ON PRACTICAL BIOLOGY

October/November 2014

The Supervisor or Teacher responsible for the subject should provide the following information.

1. Was any difficulty experienced in providing the necessary materials? If so, give brief details.

- 2. Give details of any difficulties experienced by particular candidates, giving names and candidate numbers. Reference should be made to:
 - (a) difficulties arising from faulty specimens or microscopes;
 - (b) accidents to apparatus or materials;
 - (c) assistance provided in case of colour-blindness;
 - (d) any other information that is likely to assist the Examiner, especially if this cannot be discovered from the scripts.

All other cases of individual hardship, e.g. illness or disability, should be reported direct to CIE on the normal 'Special Consideration Form' as detailed in the Handbook for Centres.

3. During the examination, the Supervisor should, **out of the sight of the candidates**, carry out **Question 1** using the same solutions and reagents as the candidates. These results should be written in the Supervisor's Report which should be enclosed with the candidates' scripts. If the scripts are in several packets, please ensure that a copy of the Supervisor's Report is enclosed with each packet of scripts.

The Invigilator should **not** carry out **Question 1**.

4. Enclose a seating plan of work benches with the scripts, giving details of the candidate numbers for the places occupied by the candidates for each session. Use separate paper for this.

Temperature of examination room°C

Results for Questions 1(a)(ii) and 1(a)(vi):

Declaration (to be signed by the Principal or the Examinations Officer)

The preparation of this practical examination has been carried out so as to maintain fully the security of the examination.

Signed Name (in block capitals) Centre number (of enclosed scripts) Centre name

If scripts are required by CIE to be despatched in more than one envelope, it is essential that a copy of the relevant Supervisor's Report and the appropriate seating plan(s) are sent inside **each envelope**.

9700/33/CI/O/N/14